



PATIENT & CAREGIVER EDUCATION

How to Set Up Your iHealth® Pulse Oximeter

This video will show you how to set up your iHealth pulse oximeter so it automatically sends your results to your MSK healthcare provider. For written instructions, read the resource *How to Set Up Your iHealth® Pulse Oximeter* (www.mskcc.org/pe/set_up_ihealth_pulse_oximeter).



Please visit

www.mskcc.org/pe/set_up_ihealth_pulse_oximeter_video
to watch this video.

For instructions for using your iHealth pulse oximeter, read the resource *How to Use Your iHealth® Pulse Oximeter* (www.mskcc.org/pe/use_ihealth_pulse_oximeter) or watch the video *How to Use Your iHealth® Pulse Oximeter* (www.mskcc.org/pe/use_ihealth_pulse_oximeter_video).

If you have any questions, contact a member of your healthcare team directly. If you're a patient at MSK and you need to reach a provider after 5:00 PM, during the weekend, or on a holiday, call 212-639-2000.

For more resources, visit www.mskcc.org/pe to search our virtual library.

How to Set Up Your iHealth® Pulse Oximeter - Last updated on October 7, 2020
All rights owned and reserved by Memorial Sloan Kettering Cancer Center