



---

PATIENT & CAREGIVER EDUCATION

# How to Disconnect Your Elastomeric Pump from Your Implanted Port

In this video, we'll show how to disconnect your elastomeric pump chemotherapy bottle from your implanted port.



Please visit [www.mskcc.org/pe/disconnect-elastomeric](http://www.mskcc.org/pe/disconnect-elastomeric) to watch this video.

For written instructions for disconnecting your elastomeric pump from your implanted port, read the “After Your Continuous Infusion” section of *Continuous Infusion with Your Elastomeric Pump* ([www.mskcc.org/pe/elastomeric-pump](http://www.mskcc.org/pe/elastomeric-pump)).

If you have questions or concerns, contact your healthcare provider. A member of your care team will answer Monday through Friday from 9 a.m. to 5 p.m. Outside those hours, you can leave a message or talk with another MSK provider. There is always a doctor or nurse on call. If you're not sure how to reach your healthcare provider, call 212-639-2000.

For more resources, visit [www.mskcc.org/pe](http://www.mskcc.org/pe) to search our virtual library.

---

How to Disconnect Your Elastomeric Pump from Your Implanted Port - Last updated on June 29, 2022

All rights owned and reserved by Memorial Sloan Kettering Cancer Center