

PATIENT & CAREGIVER EDUCATION

# How to Disconnect Your Elastomeric Pump from Your Implanted Port

In this video, we'll show how to disconnect your elastomeric pump chemotherapy bottle from your implanted port.

 Please visit [www.mskcc.org/pe/disconnect-elastomeric](http://www.mskcc.org/pe/disconnect-elastomeric) to watch this video.

For written instructions for disconnecting your elastomeric pump from your implanted port, read the “After Your Continuous Infusion” section of *Continuous Infusion with Your Elastomeric Pump* ([www.mskcc.org/pe/elastomeric-pump](http://www.mskcc.org/pe/elastomeric-pump)).

If you have questions or concerns, contact your healthcare provider. A member of your care team will answer Monday through Friday from 9 a.m. to 5 p.m. Outside those hours, you can leave a message or talk with another MSK provider. There is always a doctor or nurse on call. If you're not sure how to reach your healthcare provider, call 212-639-2000.

For more resources, visit [www.mskcc.org/pe](http://www.mskcc.org/pe) to search our virtual library.

How to Disconnect Your Elastomeric Pump from Your Implanted Port - Last updated on April 18, 2022

All rights owned and reserved by Memorial Sloan Kettering Cancer Center

## Last Updated

April 18, 2022

Learn about our [Health Information Policy](#).

## Tell us what you think

## Tell us what you think

Your feedback will help us improve the educational information we provide. Your care team cannot see anything you write on this feedback form. Please do not use it to ask about your care. If you have questions about your care, contact your healthcare provider.

While we read all feedback, we cannot answer any questions. Please do not write your name or any personal information on this feedback form.

Was this information easy to understand?

Yes

Somewhat

No

What could we have explained better?

Please do not write your name or any personal information.

Submit