



PATIENT & CAREGIVER EDUCATION

How to Set Up Your iHealth® Pulse Oximeter

This video will show you how to set up your iHealth pulse oximeter so it automatically sends your results to your MSK healthcare provider.

For written instructions, read *How To Set Up Your iHealth® Pulse Oximeter* (www.mskcc.org/pe/set_up_ihealth_pulse_oximeter).



Please visit www.mskcc.org/pe/set_up_ihealth_pulse_oximeter_video to watch this video.

For instructions for using your iHealth pulse oximeter, read *How To Use Your iHealth® Pulse Oximeter* (www.mskcc.org/pe/use_ihealth_pulse_oximeter) or watch *How to Use Your iHealth® Pulse Oximeter* (www.mskcc.org/pe/use_ihealth_pulse_oximeter_video).

If you have any questions, contact a member of your care team directly. If you're a patient at MSK and you need to reach a provider after 5 p.m., during the weekend, or on a holiday, call 212-639-2000.

For more resources, visit www.mskcc.org/pe to search our virtual library.

How to Set Up Your iHealth® Pulse Oximeter - Last updated on October 7, 2020
All rights owned and reserved by Memorial Sloan Kettering Cancer Center