



PATIENT & CAREGIVER EDUCATION

How to Set Up Your iHealth® Pulse Oximeter

This video will show you how to set up your iHealth pulse oximeter so it automatically sends your results to your MSK healthcare provider.

For written instructions, read *How To Set Up Your iHealth® Pulse Oximeter*

(www.mskcc.org/pe/set_up_ihealth_pulse_oximeter).

Please visit
 www.mskcc.org/pe/set_up_ihealth_pulse_oximeter_video
to watch this video.

For instructions for using your iHealth pulse oximeter, read *How To Use Your iHealth® Pulse Oximeter*

(www.mskcc.org/pe/use_ihealth_pulse_oximeter) or watch *How to Use Your iHealth® Pulse Oximeter*

(www.mskcc.org/pe/use_ihealth_pulse_oximeter_video).

If you have questions or concerns, contact your healthcare provider. A member of your care team will answer Monday through Friday from 9 a.m. to 5 p.m. Outside those hours, you can leave a message or talk with another MSK provider. There is always a doctor or nurse on call. If you're not sure how to reach your healthcare provider, call 212-639-2000.

For more resources, visit www.mskcc.org/pe to search our virtual library.

How to Set Up Your iHealth® Pulse Oximeter - Last updated on October 7, 2020

All rights owned and reserved by Memorial Sloan Kettering Cancer Center